

<b>GeneLab SOP for Normalizing TruSeq Stranded Total RNA Library</b>	Document No.:	GL-SOP-6.4
	Version:	1.0
	Created:	11_23_2018
	Last revised:	12_28_2018
	Last revised by:	Marie Dinh

### Purpose/Scope:

This SOP describes the steps used by NASA GeneLab to normalize TruSeq Stranded Total RNA libraries.

### Equipment:

1. TapeStation D1000

### Reagents:

1. DNase- and RNase-free water
2. TruSeq Resuspension Buffer (RSB)
3. TruSeq Total RNA library
4. Ice
5. Ice bucket
6. Low bind microtube
7. 96-well sterile plate

### Procedure:

1. Obtain average library fragment size from TapeStation D1000 without adapter dimer if average size are within similar range  $\pm 30$  bp.
2. Convert each library concentration from PicoGreen measurement ng/uL to nM using average or individual fragment size bp. (Reference PicoGreen SOP)

Equation for converting dsDNA:

$$\text{conc nM} = \frac{(\text{conc ng/uL})}{(660 \text{ g/mol} \times \text{library size bp})} \times 10^6$$

3. Dilute each library to 20nM in DNase-/RNase-free water.
4. Store at  $-20^{\circ}\text{C}$  for iSeq (GL-SOP-015.1) and NovaSeq run (GL-SOP-016.1).

Sample #	Sample	Library PicoGreen conc. (ng/ul)	Libraries (nM)	20 nM	
				Sample Vol. to 20 nM in 20ul	Water ul
1	MGS_HLU_IR_M9_D SKN_RNA_ALQ0	6.490	32.78	12.20	7.80